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**PNL-35 v5: PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY
PREPARATION USING THE ILLUMINA® DNA PREP KIT**

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
Table 2. Estimated DNA load capacity (in Mb) for Illumina reagent kits.

Organism	Estimated Genome Size (Mb)	Minimum Coverage
<i>E. coli & Shigella spp.</i>	5	40x
<i>Salmonella spp.</i>	5	30x
<i>Vibrio cholerae</i>	4	40x
<i>Vibrio spp.</i>	5	40x
<i>Listeria monocytogenes</i>	3	20x
<i>Campylobacter spp.</i>	1.6	20x
<i>Cronobacter spp.</i>	4.5	40x
<i>Yersinia spp.</i>	5	40x


Table 1. Estimated genome size (in Mb) by organism

Illumina Bead-Linked Transposome Chemistry

(A) Bead-linked transposomes mediate the simultaneous fragmentation of gDNA and the addition of Illumina sequencing primers. (B) Reduced-cycle PCR amplification amplifies sequencing ready DNA fragments and adds indexes and adapters. (C) Sequencing-ready fragments are washed and pooled.

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

Blood

Saliva

Flex Lysis Reagent Kit
Hands-on: 20 minutes
Total: 40 minutes
Reagents: BLB, PK1, SPB, EtOH

Safe Stopping Point

Flex Lysis Reagent Kit
Hands-on: 20 minutes
Total: 30 minutes
Reagents: SPB, EtOH

Safe Stopping Point

1

Tagment Genomic DNA

Hands-on: 20 minutes
Total: 40 minutes*
Reagents: BLT, TB1

2

Post Tagmentation Cleanup

Hands-on: 20 minutes
Total: 45 minutes
Reagents: TSB, TWB

3

Amplify Tagmented DNA

Hands-on: 15 minutes
Total: 45 minutes
Reagents: EPM, Index Oligos

Safe Stopping Point

4

Clean Up Libraries

Hands-on: 20 minutes
Total: 50 minutes
Reagents: RSB, SPB, EtOH

Safe Stopping Point

5

Pool Libraries

6

Whole-Genome Sequencing

On Illumina Sequencer:
iSeq, MiniSeq, MiSeq, NextSeq, HiSeq, or NovaSeq

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


*Vortex BLT mínimo 10 seg
Vortex TB1 y spin


- 1 **Tagment Genomic DNA**
Hands-on: 20 minutes
Total: 40 minutes*
Reagents: BLT, TB1
- 2 **Post Tagmentation Cleanup**
Hands-on: 20 minutes
Total: 45 minutes
Reagents: TSB, TWB
- 3 **Amplify Tagmented DNA**
Hands-on: 15 minutes
Total: 45 minutes
Reagents: EPM, Index Oligos
- 4 **Clean Up Libraries**
Hands-on: 20 minutes
Total: 50 minutes
Reagents: RSB, SPB, EtOH
- 5 **Pool Libraries**
- 6 **Whole-Genome Sequencing**
On Illumina Sequencer:
iSeq, MiniSeq, MiSeq, NextSeq, HiSeq, or NovaSeq
- 7 **Data Analysis on BaseSpace**

Safe Stopping Point


Safe Stopping Point


- 1 Tagment Genomic DNA**
Hands-on: 20 minutes
Total: 40 minutes*
Reagents: BLT, TB1
- 2 Post Tagmentation Cleanup**
Hands-on: 20 minutes
Total: 45 minutes
Reagents: TSB, TWB
- 3 Amplify Tagmented DNA**
Hands-on: 15 minutes
Total: 45 minutes
Reagents: EPM, Index Oligos
- 4 Clean Up Libraries**
Hands-on: 20 minutes
Total: 50 minutes
Reagents: RSB, SPB, EtOH
- 5 Pool Libraries**
- 6 Whole-Genome Sequencing**
On Illumina Sequencer:
iSeq, MiniSeq, MiSeq, NextSeq, HiSeq, or N
- 7 Data Analysis on BaseSpace**



Safe Stopping Point

Safe Stopping Point

LIMPIEZA BIBLIOTECAS

Organism	Instrument	Pooling volume (µl)
<i>Campylobacter</i>	MiSeq, MiniSeq	2.5
<i>Campylobacter</i>	iSeq	2.0
<i>Listeria</i>	iSeq	2.0

Table 8. Exceptions to the 5 µl pooling volume for the Illumina DNA Prep libraries

NOTE: Lower pooling volumes are required to achieve a more balanced run. This will help prevent over-coverage due to small genome size and lower coverage requirements.

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